

A 838

Key set

For setting basket heights and assembling/dismantling injector nozzles.



- Set consisting of Torx wrench and open wrench
- Suitable for Miele upper baskets and mobile units

EAN: 4002516171720 / Material number: 11054290 / Old Material Number: 69783800D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8505	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.00)
Gross weight in lbs	1 (0.00)