

A 312  
Module  
For reprocessing viscometers.



- Suitable for up to 8.0 viscometers
- Equipped with 8 x A 867 special nozzles
- Nozzle spacing 4 5/8" across the width, 4 5/8" across the depth

EAN: 4002516087397 / Material number: 11054470 / Old Material Number: 69631200D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8636 [LAB / LAB MON]	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
<b>Product category</b>	
Module for viscometers	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic, Silicone
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Viscometer [number]	8.0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	2 5/8 (66)
External dimensions, net width in inches	5 1/4 (133)
External dimensions, net depth in inches	18 7/8 (479)
External dimensions, gross height in inches	9 7/8 (250)
External dimensions, gross width in inches	9 1/16 (230)
External dimensions, gross depth in inches	19 15/16 (505)
Net weight in lbs	4 (1.58)
Gross weight in lbs	6 (2.33)