


A 836

Set of injector nozzles

For optimum reprocessing of hollow instruments used in dental applications.



- For use in A 105/1 upper basket or the A 315/1 module
- Packaging unit with 22 dental injector nozzles (A 833)

EAN: 4002516204091 / Material number: 11054660 / Old Material Number: 69783600D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
SlimLine washer-disinfectors, dental	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	4 3/4 (120)
External dimensions, gross depth in inches	6 3/4 (170)
Net weight in lbs	0 (0.20)
Gross weight in lbs	0 (0.20)