


A 313

Module

For reprocessing large-volume laboratory glassware.



- Laboratory flasks, up to 10 l
- Round flasks, up to 6 l
- Erlenmeyer flasks, up to 5 l
- Measuring flasks, up to 5 l
- Can be equipped with nozzles & holders based on individual needs
- Equipped with 5 injector positions for 2 or 3 wash items

EAN: 4002516087854 / Material number: 11055970 / Old Material Number: 69631300D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Insert/module for narrow-neck glassware	•
Module for large-volume laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	8 3/4 (221) in (mm)
External dimensions, net width in inches	9 1/4 (234) in (mm)
External dimensions, net depth in inches	19 1/2 (495) in (mm)
External dimensions, gross height in inches	10 1/4 (260)
External dimensions, gross width in inches	10 11/16 (270)
External dimensions, gross depth in inches	19 15/16 (505)
Net weight in lbs	7 (2.82)
Gross weight in lbs	8 (3.29)