

A 315 Module

For increased capacity and reprocessing of hollow instruments in dentistry.



- Increases the injector capacity from 22 to 44
- Backward-compatible with all A 151 products on the market
- Compatible with all standard dental & ENT accessories

EAN: 4002516087861 / Material number: 11055980 / Old Material Number: 69631500D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for ENT instruments	•
Insert/module for dental instruments	•
Other inserts/other modules	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Option to use a central filter	•
Application area	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Capacity	
ENT instruments [number]	22
Transmission instruments [number]	22
Siphons [number]	22
Lumened instruments [number]	22
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	145
External dimensions, net width in inches	492
External dimensions, net depth in inches	542
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	3 (1.35)
Gross weight in lbs	6 (2.35)
Accessory pack	
A 838 key set [number]	1