

A 840

Injector nozzle

With plastic supports, nozzle diameter 6, length 130 mm, 1 item.



- Suitable for mobile injector units & baskets, w/ screw connection
- Part of the SimpleLoad system
- Packaging unit with 1 nozzle
- Ideal for wash items with neck diameters between 1/2" and 3 3/8" (12 and 85 mm)
- Diameter of cap at top: 7/16" (11 mm)

EAN: 4002516188742 / Material number: 11056100 / Old Material Number: 69684000D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	5 1/8 (130) in (mm)
External dimensions, net width in inches	3 5/8 (92) in (mm)
External dimensions, net depth in inches	3 7/16 (86) in (mm)
External dimensions, gross height in inches	3 3/16 (80)
External dimensions, gross width in inches	3 3/4 (95)
External dimensions, gross depth in inches	4 3/8 (110)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.06)