

A 844

Injector nozzle

With plastic supports, nozzle diameter 2.5, length 80 mm, 1 item.



- Suitable for mobile injector units & baskets, w/ screw connection
- Ideal for wash items with neck diameters between 1/4" and 2 3/16" (6 and 55 mm)
- Diameter of cap at top: 3/16" (5 mm)
- Part of the SimpleLoad system

EAN: 4002516188780 / Material number: 11056600 / Old Material Number: 69684400D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	3 3/16 (80)
External dimensions, net width in inches	2 3/8 (60)
External dimensions, net depth in inches	2 3/16 (55)
External dimensions, gross height in inches	1 5/8 (40)
External dimensions, gross width in inches	2 9/16 (65)
External dimensions, gross depth in inches	4 11/16 (118)
Net weight in lbs	1 (0.01)
Gross weight in lbs	1 (0.03)