

A 845

Injector nozzle

With plastic supports, nozzle diameter 2.5, length 125 mm, 1 item.



- Suitable for mobile injector units & baskets, w/ screw connection
- Ideal for wash items with neck diameters between 1/4" and 2 3/16" (6 and 55 mm)
- Diameter of cap at top: 3/16" (5 mm)
- Part of the SimpleLoad system

EAN: 4002516188797 / Material number: 11056670 / Old Material Number: 69684500D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	4 15/16 (125)
External dimensions, net width in inches	2 3/8 (60)
External dimensions, net depth in inches	2 3/16 (55)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.01)
Gross weight in lbs	1 (0.01)