


A 301/4

Injector module

For optimum loading of laboratory glassware – can be equipped with 18 nozzles.

EAN: 4002516087649 / Material number: 11056780 / Old Material Number: 69630104D



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 18 injector nozzles
- Part of the SimpleLoad system
- Nozzle spacing 3" (78 mm) across the width, 3" (78 mm) across the depth
- Particularly suitable for A 842 and A 843 (nozzle diameter 4 mm)

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-neck glassware [number]	18
<b>Standard electrical connection</b>	
Number of phases	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	2 7/8 (72)
External dimensions, net width in inches	6 13/16 (173)
External dimensions, net depth in inches	18 7/8 (479)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	8 1/8 (205)
External dimensions, gross depth in inches	19 3/16 (487)
Net weight in lbs	2 (0.85)
Gross weight in lbs	5 (1.89)