


A 302/2

Injector module

For optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 32 injector nozzles
- Part of the SimpleLoad system
- Nozzle spacing 2 3/8" across the width, 2 1/4" across the depth
- Particularly suitable for A 844 and A 845 (nozzle diameter 2.5 mm)

EAN: 4002516087656 / Material number: 11056790 / Old Material Number: 69630202D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	32
Standard electrical connection	
Number of phases	0
Dimensions and weight	
External dimensions, net height in inches	2 7/8 (72)
External dimensions, net width in inches	7 11/16 (195)
External dimensions, net depth in inches	18 7/8 (479)
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	3 (1.09)
Gross weight in lbs	5 (2.16)