

A 860

Holding frames

For injector modules A 300/2 and A 300/3



- Arrangement: 2 x 4 fields
- For use on A 300/2 and A 300/3 injector modules
- For the optimum positioning of laboratory glassware
- Part of the SimpleLoad system

EAN: 4002516087663 / Material number: 11056800 / Old Material Number: 69686000D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Other inserts/other modules	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Number of compartments	8
Compartment width in inches	4 1/2 (114)
Compartment length in inches	4 1/2 (114)
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	1 3/16 (30)
External dimensions, net width in inches	9 7/8 (250)
External dimensions, net depth in inches	13 13/16 (350)
External dimensions, gross height in inches	1 7/16 (35)
External dimensions, gross width in inches	11 13/16 (300)
External dimensions, gross depth in inches	15 3/4 (400)
Net weight in lbs	1 (0.10)
Gross weight in lbs	1 (0.18)