

A 861  
Holding frames  
for A 301/4 and A 301/5 injector modules.



- Arrangement: 3 x 6 fields
- For use on A 301/4 and A 301/5 injector modules
- For the optimum positioning of laboratory glassware
- Part of the SimpleLoad system

EAN: 4002516087670 / Material number: 11056810 / Old Material Number: 69686100D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Number of compartments	18
Compartment width in inches	2 15/16 (74)
Compartment length in inches	2 15/16 (74)
Part of the SimpleLoad system	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	3 3/16 (80)
External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	13 13/16 (350)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	8 5/16 (210)
External dimensions, gross depth in inches	19 1/2 (495)
Net weight in lbs	1 (0.10)
Gross weight in lbs	5 (1.87)