

APLW 052/1

Set of injector nozzles

With plastic supports, Ø 4 mm, length 175 mm.

EAN: 4002516130918 / Material number: 11094000 / Old Material Number: 69305201D



- Suitable for e.g., lower basket APLW 095
- Packaging unit with 39 injector nozzles
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 13/64" (5 mm)

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine/accessory categorization	
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Dimensions and weight	
External dimensions, gross height in inches	7 1/2 (190)
External dimensions, gross width in inches	12 1/4 (310)
External dimensions, gross depth in inches	12 1/4 (310)
Net weight in lbs	2 (0.78)
Gross weight in lbs	3 (1.05)