


APLW 000

Mobile unit

for optimum loading of laboratory glassware on the top levels.



- Equipped with 1 integrated spray arm
- For loading inserts
- Clearance: 1' 7 1/16" x 1' 8 11/16" (484 x 525 mm)

EAN: 4002516115281 / Material number: 11104470 / Old Material Number: 69300000D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for beakers	•
Upper basket	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
Dimensions and weight	
External dimensions, net height in inches	6 5/16 (160)
External dimensions, net width in inches	22 1/16 (560)
External dimensions, net depth in inches	21 (532)
External dimensions, gross height in inches	7 1/2 (190)
External dimensions, gross width in inches	22 1/16 (560)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	9 (4.30)
Gross weight in lbs	13 (6.07)