

APLW 001

Mobile unit

for optimum loading of laboratory glassware on the bottom level.

EAN: 4002516115298 / Material number: 11104480 / Old Material Number: 69300100D



- Maximum load width: 18 1/2" (470 mm)
- Maximum load depth: 19 19/64" (490 mm)

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for beakers	•
Mobile unit for inserts	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	5 3/8 (135)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 15/16 (582)
External dimensions, gross height in inches	7 1/2 (190)
External dimensions, gross width in inches	22 1/16 (560)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	10 (4.20)
Gross weight in lbs	14 (5.90)