

APLW 003

Mobile unit

for optimum loading of large-volume laboratory glassware up to 5.2 gal. (20 l.)



**Image
in progress**

- Equipped with 4 nozzles
- For loading up to 2 bottles with a max. Ø of 11 1/32" (280 mm.)
- For loading up to 2 bottles with a max. Ø of 7 1/2" (190 mm)

EAN: 4002516115410 / Material number: 11104500 / Old Material Number: 69300300D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Flasks [number]	2
Erlenmeyer flasks [number]	2
Dimensions and weight	
External dimensions, net height in inches	400
External dimensions, net width in inches	535
External dimensions, net depth in inches	582
External dimensions, gross height in inches	22 7/8 (580)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	12 (5.08)
Gross weight in lbs	17 (7.63)