

APLW 005
 Mobile unit
 for optimum loading of up to 56 pipettes.



- For loading up to 56 pipettes arranged diagonally

EAN: 4002516115434 / Material number: 11104520 / Old Material Number: 69300500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of pipettes/vials	•
Capacity	
Pipettes [number]	56
Dimensions and weight	
External dimensions, net height in inches	12 3/16 (308)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 15/16 (582)
External dimensions, gross height in inches	22 1/2 (570)
External dimensions, gross width in inches	22 1/16 (560)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	15 (6.40)
Gross weight in lbs	21 (9.20)