

APLW 006
 Mobile unit
 for optimum loading of up to 121 pipettes.



- For loading up to 121 pipettes arranged vertically

EAN: 4002516115441 / Material number: 11104530 / Old Material Number: 69300600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of pipettes/vials	•
Capacity	
Pipettes [number]	121
Dimensions and weight	
External dimensions, net height in inches	195
External dimensions, net width in inches	535
External dimensions, net depth in inches	575
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/8 (610)
Net weight in lbs	11 (4.94)
Gross weight in lbs	17 (7.38)