

APLW 006
 Mobile unit
 for optimum loading of up to 121 pipettes.



- For loading up to 121 pipettes arranged vertically

EAN: 4002516115441 / Material number: 11104530 / Old Material Number: 69300600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of pipettes/vials	•
Capacity	
Pipettes [number]	121
Dimensions and weight	
External dimensions, net height in inches	7 11/16 (195)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 11/16 (575)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	11 (4.94)
Gross weight in lbs	17 (7.38)