

APLW 007
 Mobile unit
 for optimum loading of 2 pipette cassettes.



- For lower level
- Maximum pipette length: 20 1/2" (520 mm)

EAN: 4002516115458 / Material number: 11104540 / Old Material Number: 69300700D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of pipettes/vials	•
Dimensions and weight	
External dimensions, net height in inches	22 7/8 (580)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	16 9/16 (420)
External dimensions, gross height in inches	22 7/8 (580)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	16 (7.18)
Gross weight in lbs	23 (9.98)