

APLW 007
 Mobile unit
 for optimum loading of 2 pipette cassettes.



- For lower level
- Maximum pipette length: 20 1/2" (520 mm)

EAN: 4002516115458 / Material number: 11104540 / Old Material Number: 69300700D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of pipettes/vials	•
Dimensions and weight	
External dimensions, net height in inches	580
External dimensions, net width in inches	535
External dimensions, net depth in inches	420
External dimensions, gross height in inches	22 7/8 (580)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/8 (610)
Net weight in lbs	16 (7.18)
Gross weight in lbs	23 (9.98)