


APLW 008  
 Mobile unit  
 for optimum loading of 3 pipette cassettes.



- For lower level
- Pipette cassette dimensions: W 3 1/4" x H 15 1/8" x D 3 1/4"

EAN: 4002516115465 / Material number: 11104550 / Old Material Number: 69300800D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6111	•
<b>Product category</b>	
Mobile unit for pipettes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of pipettes/vials	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	18 (456)
External dimensions, net width in inches	7 1/16 (179)
External dimensions, net depth in inches	14 7/16 (366)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	17 (7.64)
Gross weight in lbs	22 (9.64)