


APLW 008  
 Mobile unit  
 for optimum loading of 3 pipette cassettes.



- For lower level
- Pipette cassette dimensions: W 3 1/4" x H 15 1/8" x D 3 1/4"

EAN: 4002516115465 / Material number: 11104550 / Old Material Number: 69300800D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6111	•
<b>Product category</b>	
Mobile unit for pipettes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of pipettes/vials	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	456
External dimensions, net width in inches	179
External dimensions, net depth in inches	366
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/8 (610)
Net weight in lbs	17 (7.64)
Gross weight in lbs	22 (9.64)