

APLW 012
Upper basket
with 84 empty injector spaces for nozzles Ø 3/32" (2.5 mm.)



- Only suitable for use in the upper levels

EAN: 4002516115502 / Material number: 11104590 / Old Material Number: 69301200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	84
Dimensions and weight	
External dimensions, net height in inches	6 1/4 (158)
External dimensions, net width in inches	21 (532)
External dimensions, net depth in inches	22 1/2 (570)
External dimensions, gross height in inches	13 13/16 (350)
External dimensions, gross width in inches	22 1/16 (560)
External dimensions, gross depth in inches	23 5/8 (600)
Net weight in lbs	11 (4.60)
Gross weight in lbs	17 (7.40)