

APLW 013/2

Upper basket

equipped with 32 APLW 051 injector nozzles and 10 APLW 055 injector nozzles



**Image
in progress**

- Only suitable for use in the upper levels

EAN: 4002516115533 / Material number: 11104620 / Old Material Number: 69301302D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in inches	9 7/8 (250)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 1/2 (570)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	11 (4.55)
Gross weight in lbs	20 (9.00)
Accessory pack	
APLW 051, injector nozzle Ø 3/16" (4 mm), length 4 5/16" (110 mm) [number]	32
APLW 055, injector nozzle Ø 1/4" (6 mm), length 6 7/8" (175 mm) [number]	10