

APLW 014  
Upper basket  
with 20 empty injector spaces.



- Only suitable for use in the upper levels

EAN: 4002516115540 / Material number: 11104630 / Old Material Number: 69301400D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-neck glassware [number]	20
Flasks [number]	20
Erlenmeyer flasks [number]	20
Round flasks [number]	20
<b>Dimensions and weight</b>	
External dimensions, net height in inches	0 7/8 (21)
External dimensions, net width in inches	0 13/16 (20)
External dimensions, net depth in inches	0 11/16 (16)
External dimensions, gross height in inches	13 7/16 (340)
External dimensions, gross width in inches	21 7/8 (555)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	10 (4.22)
Gross weight in lbs	15 (6.65)