

APLW 014/1  
Upper basket  
equipped with 20 APLW 055 injector nozzles



- Only suitable for use in the upper levels

EAN: 4002516115557 / Material number: 11104640 / Old Material Number: 69301401D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
<b>Capacity</b>	
Narrow-neck glassware [number]	20
Flasks [number]	20
Erlenmeyer flasks [number]	20
Round flasks [number]	20
<b>Dimensions and weight</b>	
External dimensions, net height in inches	9 7/8 (250)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 1/2 (570)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	9 (4.22)
Gross weight in lbs	19 (8.46)
<b>Accessory pack</b>	
APLW 055, injector nozzle Ø 1/4" (6 mm), length 6 7/8" (175 mm) [number]	20