

APLW 015
Upper basket
with 16 empty injector spaces.



- Only suitable for use in the upper levels

EAN: 4002516115571 / Material number: 11104650 / Old Material Number: 69301500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	16
Flasks [number]	16
Erlenmeyer flasks [number]	16
Round flasks [number]	16
Dimensions and weight	
External dimensions, net height in inches	160
External dimensions, net width in inches	532
External dimensions, net depth in inches	560
External dimensions, gross height in inches	13 1/2 (340)
External dimensions, gross width in inches	22 1/8 (560)
External dimensions, gross depth in inches	23 5/8 (600)
Net weight in lbs	9 (4.00)
Gross weight in lbs	15 (6.60)