

APLW 018  
Upper basket  
with 121 + 24 empty injector spaces.



- Only suitable for use in the upper levels

EAN: 4002516115601 / Material number: 11104680 / Old Material Number: 69301800D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
Narrow-neck glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in inches	6 7/16 (162)
External dimensions, net width in inches	21 1/16 (534)
External dimensions, net depth in inches	22 7/16 (569)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	12 (5.33)
Gross weight in lbs	18 (7.76)