


APLW 018/1

Upper basket

equipped with 121 APLW 048 injector nozzles and 24 APLW 051 injector nozzles.

EAN: 4002516115717 / Material number: 11104690 / Old Material Number: 69301801D



- Only suitable for use in the upper levels

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
Narrow-neck glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in inches	7 7/8 (200)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 1/2 (570)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	16 (6.92)
Gross weight in lbs	26 (11.42)
<b>Accessory pack</b>	
APLW 051, injector nozzle Ø 3/16" (4 mm), length 4 5/16" (110 mm) [number]	32
APLW 048 injector nozzle Ø 1/8" (2.5 mm), length 3 1/8" (80 mm) [number]	121