

APLW 019
 Upper basket
 with 121 empty injector spaces for nozzles Ø 3/32" (2.5 mm.)



**Image
 in progress**

- Only suitable for use in the upper levels

EAN: 4002516115724 / Material number: 11104700 / Old Material Number: 69301900D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Phials [number]	121
Dimensions and weight	
External dimensions, net height in inches	93
External dimensions, net width in inches	534
External dimensions, net depth in inches	569
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	13 (5.52)
Gross weight in lbs	18 (7.96)