

APLW 019

Upper basket

with 121 empty injector spaces for nozzles Ø 3/32" (2.5 mm.)



Image
in progress

- Only suitable for use in the upper levels

EAN: 4002516115724 / Material number: 11104700 / Old Material Number: 69301900D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Phials [number]	121
Dimensions and weight	
External dimensions, net height in inches	3 11/16 (93)
External dimensions, net width in inches	21 1/16 (534)
External dimensions, net depth in inches	22 7/16 (569)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	13 (5.52)
Gross weight in lbs	18 (7.96)