


APLW 019/1  
Upper basket  
equipped with 121 APLW 047 injector nozzles.



- Only suitable for use in the upper levels

EAN: 4002516115731 / Material number: 11104710 / Old Material Number: 69301901D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
<b>Dimensions and weight</b>	
External dimensions, net height in inches	200
External dimensions, net width in inches	535
External dimensions, net depth in inches	570
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	15 (6.64)
Gross weight in lbs	20 (9.06)
<b>Companion pack</b>	
APLW 047 injector nozzle Ø 1/8" (2.5 mm), length 2" (50 mm) [number]	121