

APLW 022  
 Lower basket  
 with 70 empty injector spaces for nozzles Ø 3/32" (2.5 mm.)



- Only suitable for use in the lowest level

EAN: 4002516115779 / Material number: 11104750 / Old Material Number: 69302200D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-neck glassware [number]	70
Flasks [number]	70
Erlenmeyer flasks [number]	70
Round flasks [number]	70
<b>Dimensions and weight</b>	
External dimensions, net height in inches	580
External dimensions, net width in inches	515
External dimensions, net depth in inches	534
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (562)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	13 (5.83)
Gross weight in lbs	19 (8.27)