

APLW 023  
 Lower basket  
 with 56 empty injector spaces for nozzles Ø 3/32" (2.5 mm.)



- Only suitable for use in the lowest level

EAN: 4002516115786 / Material number: 11104760 / Old Material Number: 69302300D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-neck glassware [number]	56
Flasks [number]	56
Erlenmeyer flasks [number]	56
Round flasks [number]	56
<b>Dimensions and weight</b>	
External dimensions, net height in inches	22 11/16 (575)
External dimensions, net width in inches	19 7/8 (504)
External dimensions, net depth in inches	18 15/16 (480)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 3/16 (562)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	13 (5.48)
Gross weight in lbs	18 (7.91)