

APLW 024  
Lower basket  
with 42 empty injector spaces for nozzles.



- Only suitable for use in the lowest level

EAN: 4002516115793 / Material number: 11104770 / Old Material Number: 69302400D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in inches	20
External dimensions, net width in inches	21
External dimensions, net depth in inches	22
External dimensions, gross height in inches	13 7/8 (352)
External dimensions, gross width in inches	22 1/4 (562)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	12 (5.22)
Gross weight in lbs	17 (7.66)