

APLW 024/2

Lower basket

equipped with 32 APLW 051 injector nozzles and 10 APLW 055 injector nozzles.



**Image  
in progress**

- Only suitable for use in the lowest level

EAN: 4002516115816 / Material number: 11104790 / Old Material Number: 69302402D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, Black, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-neck glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
<b>Dimensions and weight</b>	
External dimensions, net height in inches	8 11/16 (220)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 11/16 (575)
External dimensions, gross height in inches	13 7/8 (352)
External dimensions, gross width in inches	22 3/16 (562)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	14 (5.97)
Gross weight in lbs	24 (10.48)
<b>Accessory pack</b>	
APLW 051, injector nozzle Ø 3/16" (4 mm), length 4 5/16" (110 mm) [number]	32
APLW 055, injector nozzle Ø 1/4" (6 mm), length 6 7/8" (175 mm) [number]	10