

APLW 024/3

Lower basket

equipped with 30 APLW 066 injector nozzles and 12 APLW 067 injector nozzles.



**Image
in progress**

- Only suitable for use in the lowest level

EAN: 4002516115823 / Material number: 11104800 / Old Material Number: 69302403D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in inches	8 11/16 (220)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 11/16 (575)
External dimensions, gross height in inches	13 7/8 (352)
External dimensions, gross width in inches	22 3/16 (562)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	12 (5.22)
Gross weight in lbs	31 (13.87)
Accessory pack	
APLW 066 injector nozzle Ø 1/4" (6 mm), length 6 7/8" (175 mm) [number]	30
APLW 067 injector nozzle Ø 1/4" (6 mm), length 8 7/8" (225 mm) [number]	12