

APLW 025/1

Lower basket

equipped with 12 APLW 056 injector nozzles and 12 APLW 052 injector nozzles.



**Image
in progress**

- Only suitable for use in the lowest level
- 9 1/16" x 19 19/64" (230 x 490 mm.) usable space

EAN: 4002516115847 / Material number: 11104820 / Old Material Number: 69302501D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
Dimensions and weight	
External dimensions, net height in inches	10 11/16 (270)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 1/2 (570)
External dimensions, gross height in inches	13 13/16 (350)
External dimensions, gross width in inches	22 1/16 (560)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	12 (5.25)
Gross weight in lbs	19 (8.39)
Accessory pack	
APLW 052 injector nozzle Ø 3/16" (4 mm), length 6 7/8" (175 mm) [number]	12
APLW 056 injector nozzle Ø 1/4" (6 mm), length 8 7/8" (225 mm) [number]	12