

APLW 025/1

Lower basket

equipped with 12 APLW 056 injector nozzles and 12 APLW 052 injector nozzles.



**Image
in progress**

- Only suitable for use in the lowest level
- 9 1/16" x 19 19/64" (230 x 490 mm.) usable space

EAN: 4002516115847 / Material number: 11104820 / Old Material Number: 69302501D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
Dimensions and weight	
External dimensions, net height in inches	270
External dimensions, net width in inches	535
External dimensions, net depth in inches	570
External dimensions, gross height in inches	13 7/8 (350)
External dimensions, gross width in inches	22 1/8 (560)
External dimensions, gross depth in inches	24 1/8 (610)
Net weight in lbs	12 (5.25)
Gross weight in lbs	19 (8.39)
Companion pack	
APLW 052 injector nozzle Ø 3/16" (4 mm), length 6 7/8" (175 mm) [number]	12
APLW 056 injector nozzle Ø 1/4" (6 mm), length 8 7/8" (225 mm) [number]	12