


APLW 026/1
Lower basket
equipped with 20 APLW 056 injector nozzles.

EAN: 4002516115861 / Material number: 11104840 / Old Material Number: 69302601D



- Only suitable for use in the lowest level

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	20
Flasks [number]	20
Erlenmeyer flasks [number]	20
Round flasks [number]	20
Dimensions and weight	
External dimensions, net height in inches	270
External dimensions, net width in inches	535
External dimensions, net depth in inches	570
External dimensions, gross height in inches	13 1/2 (340)
External dimensions, gross width in inches	21 7/8 (555)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	11 (4.92)
Gross weight in lbs	21 (9.20)
Companion pack	
APLW 056 injector nozzle Ø 1/4" (6 mm), length 8 7/8" (225 mm) [number]	20