


APLW 026/1
Lower basket
equipped with 20 APLW 056 injector nozzles.

EAN: 4002516115861 / Material number: 11104840 / Old Material Number: 69302601D



- Only suitable for use in the lowest level

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	20
Flasks [number]	20
Erlenmeyer flasks [number]	20
Round flasks [number]	20
Dimensions and weight	
External dimensions, net height in inches	10 11/16 (270)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 1/2 (570)
External dimensions, gross height in inches	13 7/16 (340)
External dimensions, gross width in inches	21 7/8 (555)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	11 (4.92)
Gross weight in lbs	21 (9.20)
Accessory pack	
APLW 056 injector nozzle Ø 1/4" (6 mm), length 8 7/8" (225 mm) [number]	20