

APLW 029
 Lower basket
 with 27 empty injector spaces for nozzles.



- Only suitable for use in the lowest level

EAN: 4002516115892 / Material number: 11104870 / Old Material Number: 69302900D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	12
Flasks [number]	12
Erlenmeyer flasks [number]	12
Round flasks [number]	12
Dimensions and weight	
External dimensions, net height in inches	5 7/16 (137)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 15/16 (582)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	11 (4.74)
Gross weight in lbs	16 (7.18)