

APLW 029
 Lower basket
 with 27 empty injector spaces for nozzles.



- Only suitable for use in the lowest level

EAN: 4002516115892 / Material number: 11104870 / Old Material Number: 69302900D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	12
Flasks [number]	12
Erlenmeyer flasks [number]	12
Round flasks [number]	12
Dimensions and weight	
External dimensions, net height in inches	137
External dimensions, net width in inches	535
External dimensions, net depth in inches	582
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	11 (4.74)
Gross weight in lbs	16 (7.18)