

APLW 030
 Lower basket
 with 12 empty injector spaces for nozzles.



- Only suitable for use in the lowest level
- 8 21/32" x 19 19/64" (220 x 490 mm.) usable space

EAN: 4002516115908 / Material number: 11104880 / Old Material Number: 69303000D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	12
Flasks [number]	12
Erlenmeyer flasks [number]	12
Round flasks [number]	12
Dimensions and weight	
External dimensions, net height in inches	5 3/8 (135)
External dimensions, net width in inches	16 3/4 (424)
External dimensions, net depth in inches	22 15/16 (582)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	10 (4.32)
Gross weight in lbs	15 (6.72)