

APLW 031  
 Lower basket  
 with 9 empty injector spaces for nozzles.



- Only suitable for use in the lowest level

EAN: 4002516115922 / Material number: 11104900 / Old Material Number: 69303100D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-neck glassware [number]	9
Flasks [number]	9
Erlenmeyer flasks [number]	9
Round flasks [number]	9
<b>Dimensions and weight</b>	
External dimensions, net height in inches	5 7/16 (137)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	22 15/16 (582)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/4 (565)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	10 (4.52)
Gross weight in lbs	16 (6.95)