

APLW 032  
 Lower basket  
 with 121 + 24 empty injector spaces for nozzles.



**Image  
 in progress**

- Only suitable for use in the lowest level

EAN: 4002516115939 / Material number: 11104910 / Old Material Number: 69303200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Phials [number]	121
Narrow-neck glassware [number]	24
Flasks [number]	24
Dimensions and weight	
External dimensions, net height in inches	22 1/2 (570) in (mm)
External dimensions, net width in inches	21 1/8 (536) in (mm)
External dimensions, net depth in inches	2 13/16 (70) in (mm)
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/16 (560)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	14 (5.99)
Gross weight in lbs	19 (8.31)