

APLW 032/1

Lower basket

equipped with 121 APLW 048 injector nozzles and 24 APLW 056 injector nozzles.



**Image**  
in progress

- Only suitable for use in the lowest level

EAN: 4002516115946 / Material number: 11104920 / Old Material Number: 69303201D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Phials [number]	121
Narrow-neck glassware [number]	24
Flasks [number]	24
Dimensions and weight	
External dimensions, net height in inches	270
External dimensions, net width in inches	535
External dimensions, net depth in inches	575
External dimensions, gross height in inches	14 (355)
External dimensions, gross width in inches	22 1/8 (560)
External dimensions, gross depth in inches	24 1/8 (610)
Net weight in lbs	14 (5.99)
Gross weight in lbs	31 (13.98)
Companion pack	
APLW 048 injector nozzle Ø 1/8" (2.5 mm), length 3 1/8" (80 mm) [number]	121
APLW 056 injector nozzle Ø 1/4" (6 mm), length 8 7/8" (225 mm) [number]	24