


APLW 034

Lid

height 9 27/32" (250 mm), for fixing delicate wash items.



EAN: 4002516115960 / Material number: 11104940 / Old Material Number: 69303400D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other inserts	•
Cover	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	9 7/8 (250) in (mm)
External dimensions, net width in inches	9 13/16 (248) in (mm)
External dimensions, net depth in inches	9 13/16 (248) in (mm)
External dimensions, gross height in inches	13 (330)
External dimensions, gross width in inches	12 5/16 (312)
External dimensions, gross depth in inches	12 7/16 (315)
Net weight in lbs	2 (0.52)
Gross weight in lbs	4 (1.38)