

APLW 036

Lid

height 1 49/64" (45 mm), for fixing delicate wash items.



EAN: 4002516115984 / Material number: 11104960 / Old Material Number: 69303600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other inserts	•
Cover	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	248
External dimensions, net width in inches	248
External dimensions, net depth in inches	45
External dimensions, gross height in inches	11 (278)
External dimensions, gross width in inches	2 1/4 (57)
External dimensions, gross depth in inches	11 1/4 (285)
Net weight in lbs	1 (0.39)
Gross weight in lbs	2 (0.53)